

## **Experimental Study of a Microfocus X-ray Generator Coupled to a Confocal Max-Flux<sup>TM</sup> Optic for Protein Crystallography**

Boris Verman, Bonglea Kim, Licai Jiang,  
John McGill, Nick Grupido, Doug Wilcox, Rick Smith  
Osmic Inc., 1788 Northwood Drive, Troy, MI 48084

Cheng Yang, Adam Courville, Charles N. Stence, Joseph D. Ferrara,  
Rigaku/Molecular Structure Corporation,  
9009 New Trails Drive, The Woodlands, TX, 77381-5209

Jimpei Harada, Masaru Kuribayashi, Kazuhiko Omote, Akihito Yamano  
Rigaku Corporation,  
3-9-12 Matsubara-Cho, Akishima-Shi, Tokyo 196-8666, Japan

In an x-ray diffraction measurement system, the combination of x-ray source and optics determines the beam characteristics. Particularly, the capture angle of the optic and the source usage efficiency define the total flux. A gradient d-spacing multilayer, with its configurable d-spacing distribution and generally larger reflection angle than a total reflection mirror, is an excellent choice for obtaining a large capture angle. On the other hand, high source usage efficiency with a large source, while yielding higher flux, also increases divergence and spectral background. Thus, high source usage efficiency may not be desirable for some applications. A high-brilliance microfocus x-ray source combined with a specially designed multilayer optic should be an ideal device, offering the following performance characteristics: high flux, low divergence, and low spectral background.


With geometric modeling and ray tracing simulations<sup>1</sup> we have designed, fabricated and tested a prototype system. The system includes a microfocus x-ray generator manufactured by Bede Scientific and a two-dimensional focusing optic manufactured by Osmic, the Microfocus Confocal Max-Flux<sup>TM</sup> Optic ( $\mu$ CMF). The combined system is called the MicroMax<sup>TM</sup>. A schematic drawing of the system is shown in Figure 1. The system parameters are given in Table 1.

In this study we compare the performance of the MicroMax x-ray source with the

---

<sup>1</sup> Configuration study of Confocal Max-Flux<sup>TM</sup> Optical System by Using a Ray Tracing Method. Licai Jiang, Boris Verman, Karsten Dan Joensen. P095 at this conference.

current benchmark for the home laboratory: the Blue-3 system consisting of a Rigaku RU-3HR with an Osmic CMF12-38Cu6 optic. This configuration has been described elsewhere<sup>2</sup>.




**Figure 1.** Schematic drawing of a focusing x-ray system

**Table1.** System parameters of the prototype system.

| System parameters                 |      |
|-----------------------------------|------|
| Source size ( $\mu\text{m}$ )     | 20   |
| Source power (W)                  | 24   |
| Source-optic distance (mm)        | 65   |
| Source-focus distance (mm)        | 700  |
| Length of the optic (mm)          | 80   |
| Capture angle ( $^\circ$ )        | 2.06 |
| Convergent angle (mR)             | 2.91 |
| Center d-spacing ( $\text{\AA}$ ) | 35   |


<sup>2</sup> Yang, C., Courville, A., Ferrara, J. (1999) *Acta Crystallographica*, **D55**, 1681-1689.



**Figure 2.** MicroMax system consisting of a Bede Microsource x-ray generator and Osmic  $\mu$ CMF optic.

### **Experimental:**

The  $\mu$ CMF optic designed for a microfocus source has more stringent engineering requirements as compared to the requirements for a CMF optic used with a rotating anode generator. The gradient of the multilayer d-spacing and the gradient of the local radius of curvature of the multilayer are much larger, while the average radius is much smaller than that of the optic used with a rotating anode. Figure 3 shows the comparison between the Confocal MaxFlux (CMF) optic designed for a rotating anode, specifically the Blue optic, and a CMF optic designed specifically for a microfocus x-ray generator,  $\mu$ CMF or mCMF.



**Figure 3.** Comparison between the Blue optic and the iCMF optic.

In order to assess the properties of the MicroMax and to compare those properties to the Blue-3 system we measured the spectral purity, beam profile at the crystal position, divergence, useable flux and diffraction on the same lysozyme crystal. We used two configurations of the MicroMax. In the first configuration we set limiting apertures of 0.3 mm, 0.5 mm and 1.0 mm diameter at a distance of 495 mm from the source. The sample position was 505 mm from the source. The purpose of this experiment was to see how the MicroMax performed in a configuration optimized for the measurement of samples with long unit cells.

In the second configuration, the high brilliancy configuration, a 0.5 mm diameter limiting aperture was set 690 mm from the source and the sample was placed 700 mm from the source. This configuration puts the sample near the focal plane of the optic, providing the maximum flux on a small sample.

## Results:

The results of the measurements of the physical properties of the beam are provided in Table 2, 3 and 4 and Figures 4-8. The results of data collection on a single frozen lysozyme crystal of dimensions 0.40 mm x 0.25 mm x 0.10 mm collected with 0.5 degree oscillations of 2 minutes on an R-Axis IV++ detector are shown in Tables 5-8. Lastly Figure 9 shows the resolution of a 421 Å axis with

the MicroMax at 450 mm crystal-to-detector distance and Table 9 provides the processing results.

### **Discussion:**

As with all optics systems based on multilayer technology the spectral purity is quite good for both the Blue-3 and MicroMax systems, typically 98% or better. Spectral purity for total reflection systems is typically much worse, with 10% or more, white radiation.

The useable flux for the long axis configuration (LAC) is about 80% of the high brilliancy configuration (HBC). The results of the processing of the lysozyme data suggest that the useable flux as seen by the crystal is 18% for the LAC as compared to the HBC. This requires reinvestigation since the difference of 18% and 80% is greater than one would expect from the error in the experiment. The predicted theoretical flux<sup>1</sup> for a 0.5 mm aperture for the LAC and HBC is  $1.5 \times 10^8$  and  $2.3 \times 10^8$  photons per second, which compares favorably with the observed intensities of  $2.0 \times 10^8$  and  $2.2 \times 10^8$  photons per second.

The Blue-3 system provides about twice the useable flux for all sample sizes as compared to the MicroMax. Data collection on the lysozyme crystal with the HBC of the MicroMax compares favorably to data collected with the Blue-3 system. The average intensity of the MicroMax data set is 46% of the Blue-3 data set, which is consistent with the useable flux measurements for 0.3 and 0.5 mm apertures. The  $R_{\text{merge}}$  for the Blue-3 is 0.036 and for MicroMax HBC data it is 0.040. The expected  $R_{\text{merge}}$  for the HBC data would be about 0.050 based on the change in counting statistics. The reflection size is not likely to be a factor since it is nearly the same for both data sets, see Figures 4 and 8. Clearly, effects other than counting statistics alone are present and this requires further investigation.

Figure 9 shows that in the LAC, resolution of a 421 Å axis is possible; the reflections are resolved well enough to process the data as shown in Table 9. This was not possible with the MicroMax in the high brilliancy configuration or for the Blue-3 system.

## Conclusion:

In the case of a routine data collection the MicroMax performs at about 46% the level of the laboratory reference Blue-3 system using useable flux as the metric of performance. Using data quality as the metric the MicroMax performs at 90%. In the case of a difficult sample with a long unit cell the MicroMax could be configured to allow accurate data collection; the Blue-3 system could not.

Future developments will include the use of a 50 W source running at 40 W, providing a 60% increase in power over the 24 W loading used for this prototype. The increase in flux should scale linearly with the loading; the performance of the MicroMax should be nearly equivalent to or surpass that of the Blue-3 system.

**Table 2.** Comparison of beam properties for the Blue-3 system and the MicroMax system with a source-to-sample distance of 505 mm and limiting apertures of 0.3 mm, 0.5 mm and 1.0 mm at 495 mm source-to-sample distance and a source-to-sample distance of 700 mm and a 0.5 mm limiting aperture at 690 mm.


| System                                   | Blue-3<br>System<br>with a 1.0 x<br>0.5 mm<br>collimator | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.3 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>1.0 mm<br>aperture | MicroMax<br>System with<br>700 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture |
|--|--|---|---|---|---|
| Properties                               |  |   |   |   |   |
| Horizontal<br>FWHM (mm)                  | 0.46(1)  | 0.190(9)  | 0.29(1)   | 0.358(4)  | 0.246(5)  |
| Vertical<br>FWHM (mm)                    | 0.38(1)  | 0.20(1)   | 0.31(2)   | 0.405(7)  | 0.31(1)   |
| Mean<br>Horizontal<br>Divergence<br>(mR) | 2.3(1)   | 1.19(5)   | 0.84(5)   | 1.2(1)  | 2.28(4)   |
| Mean<br>Vertical<br>Divergence<br>(mR)   | 2.7(1)   | 2.09(4)   | 1.45(5)   | 1.24(6)   | 2.32(5)   |
| CuK <sub>α</sub> (%)                     | 97.7   | 98.3  | 97.8  | 97.00   | -----   |
| FeK <sub>α</sub> (%)                     | 0.22   | 0.21  | 0.25  | 0.22  | -----   |
| White<br>radiation (%)                   | 2.10   | 1.46  | 1.91  | 2.78  | -----   |

**Table 3.** Comparison of useable flux in pin diode units for the Blue-3 system and the MicroMax system with a source-to-sample distance of 505 mm and limiting apertures of 0.3 mm, 0.5 mm and 1.0 mm at 495 mm source-to-sample distance and a source-to-sample distance of 700 mm and a 0.5 mm limiting aperture at 690 mm.


| System<br><br>Aperture<br>(mm) | Blue-3<br>System<br>with a 1.0 x<br>0.5 mm<br>collimator | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.3 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>1.0 mm<br>aperture | MicroMax<br>System with<br>700 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture |
|--------------------------------|--|---|---|---|---|
| 1.2                            | 10.51  | 1.75  | 3.96  | 6.11  | 4.57  |
| 1.0                            | 10.52  | 1.75  | 3.96  | 6.11  | 4.57  |
| 0.6                            | 10.41  | 1.75  | 3.96  | 5.12  | 4.57  |
| 0.5                            | 10.22  | 1.75  | 3.96  | 4.39  | 4.34  |
| 0.3                            | 4.51   | 1.60  | 1.66  | 1.82  | 2.31  |
| 0.2                            | 1.69   | 0.89  | 0.91  | 0.91  | 1.12  |
| 0.1                            | 0.45   | 0.16  | 0.22  | 0.19  | 0.18  |

**Table 4.** Comparison of useable flux in  $10^6$  photons/second for the Blue-3 system and the MicroMax system with a source-to-sample distance of 505 mm and limiting apertures of 0.3 mm, 0.5 mm and 1.0 mm at 495 mm source-to-sample distance and a source-to-sample distance of 700 mm and a 0.5 mm limiting aperture at 690 mm.

| System<br><br>Aperture<br>(mm) | Blue-3<br>System<br>with a 1.0 x<br>0.5 mm<br>collimator | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.3 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture | MicroMax<br>System with<br>505 mm<br>source-to-<br>sample<br>distance and<br>1.0 mm<br>aperture | MicroMax<br>System with<br>700 mm<br>source-to-<br>sample<br>distance and<br>0.5 mm<br>aperture |
|--------------------------------|--|---|---|---|---|
| 1.2                            | 537  | 90  | 202   | 313   | 234   |
| 1.0                            | 538  | 90  | 202   | 313   | 234   |
| 0.6                            | 532  | 90  | 202   | 262   | 234   |
| 0.5                            | 523  | 90  | 202   | 225   | 222   |
| 0.3                            | 230  | 82  | 85  | 93  | 118   |
| 0.2                            | 86   | 46  | 47  | 47  | 57  |
| 0.1                            | 23   | 8.2   | 11.2  | 9.7   | 9.8   |




**Figure 4.** Divergence plots in the horizontal and vertical directions for the Blue-3 system.




**Figure 5.** Divergence plots in the horizontal and vertical directions for the MicroMax with a source-to-sample distance of 505 mm and a 0.3 mm aperture.






**Figure 6.** Divergence plots in the horizontal and vertical directions for the MicroMax with a source-to-sample distance of 505 mm and a 0.5 mm aperture.



**Figure 7.** Divergence plots in the horizontal and vertical directions for the MicroMax with a source-to-sample distance of 505 mm and a 1.0 mm aperture.



**Figure 8.** Divergence plots in the horizontal and vertical directions for the MicroMax with a source-to-sample distance of 700 mm and a 0.5 mm aperture.

**Table 6.** Processing results for lysozyme collected with the Blue-3 system.


| Resolution<br>range | Average<br>counts | Num<br>obs | Num<br>rejs | Num<br>ovlps | Num<br>mults | <I/<br>sig> | ChiSq<br>norm | Rmerge<br>shell | Rmerge<br>cumul |
|---------------------|-------------------|------------|-------------|--------------|--------------|-------------|---------------|-----------------|-----------------|
| 24.93 - 4.09        | 102786            | 7353       | 3           | 7343         | 1036         | 31.7        | 0.41          | 0.025           | 0.025           |
| 4.09 - 3.25         | 124492            | 7417       | 3           | 7402         | 949          | 33.2        | 0.37          | 0.026           | 0.026           |
| 3.25 - 2.84         | 68252             | 7442       | 5           | 7424         | 924          | 32.6        | 0.56          | 0.032           | 0.027           |
| 2.84 - 2.58         | 42059             | 7365       | 13          | 7338         | 899          | 32.1        | 0.69          | 0.036           | 0.028           |
| 2.58 - 2.39         | 36094             | 7320       | 15          | 7288         | 885          | 31.9        | 0.76          | 0.038           | 0.029           |
| 2.39 - 2.25         | 31711             | 7255       | 32          | 7200         | 868          | 31.5        | 1.09          | 0.048           | 0.031           |
| 2.25 - 2.14         | 26267             | 7212       | 27          | 7169         | 870          | 30.9        | 1.60          | 0.066           | 0.033           |
| 2.14 - 2.05         | 22020             | 7014       | 36          | 6958         | 845          | 30.0        | 1.92          | 0.074           | 0.035           |
| 2.05 - 1.97         | 16905             | 4666       | 58          | 4541         | 739          | 23.9        | 2.24          | 0.080           | 0.036           |
| 1.97 - 1.90         | 13810             | 1613       | 30          | 1495         | 383          | 17.9        | 2.57          | 0.087           | 0.036           |
| 24.93 - 1.90        | 53335             | 64657      | 222         | 64158        | 8398         | 30.3        | 1.04          | 0.036           | 0.036           |

**Table 7.** Processing results for lysozyme collected with the MicroMax in the long axis configuration with a 0.5 mm limiting aperture.

| Resolution<br>range | Average<br>counts | Num<br>obs | Num<br>rejs | Num<br>ovlps | Num<br>mults | <I/<br>sig> | ChiSq<br>norm | Rmerge<br>shell | Rmerge<br>cumul |
|---------------------|-------------------|------------|-------------|--------------|--------------|-------------|---------------|-----------------|-----------------|
| 24.96 - 4.09        | 8484              | 7347       | 8           | 7332         | 1036         | 28.5        | 0.57          | 0.033           | 0.033           |
| 4.09 - 3.25         | 10230             | 7465       | 11          | 7442         | 955          | 29.4        | 0.71          | 0.038           | 0.036           |
| 3.25 - 2.84         | 5603              | 7388       | 12          | 7363         | 920          | 26.7        | 0.84          | 0.048           | 0.038           |
| 2.84 - 2.58         | 3379              | 7343       | 12          | 7317         | 901          | 24.0        | 0.94          | 0.061           | 0.041           |
| 2.58 - 2.39         | 2922              | 7317       | 16          | 7284         | 886          | 23.0        | 1.02          | 0.071           | 0.044           |
| 2.39 - 2.25         | 2537              | 7237       | 36          | 7178         | 868          | 21.7        | 1.27          | 0.090           | 0.047           |
| 2.25 - 2.14         | 2068              | 7296       | 36          | 7244         | 875          | 20.0        | 1.47          | 0.110           | 0.051           |
| 2.14 - 2.05         | 1696              | 7024       | 46          | 6958         | 847          | 18.0        | 1.50          | 0.127           | 0.054           |
| 2.05 - 1.97         | 1263              | 4578       | 32          | 4477         | 733          | 12.9        | 1.54          | 0.143           | 0.056           |
| 1.97 - 1.90         | 976               | 1570       | 19          | 1463         | 373          | 8.5         | 1.38          | 0.154           | 0.057           |
| 24.96 - 1.90        | 4335              | 64565      | 228         | 64058        | 8394         | 22.2        | 1.08          | 0.057           | 0.057           |

**Table 8.** Processing results for lysozyme collected with the MicroMax in the high brilliancy configuration.

| Resolution<br>range | Average<br>counts | Num<br>obs | Num<br>rejs | Num<br>ovlps | Num<br>mults | <I/<br>sig> | ChiSq<br>norm | Rmerge<br>shell | Rmerge<br>cumul |
|---------------------|-------------------|------------|-------------|--------------|--------------|-------------|---------------|-----------------|-----------------|
| 24.91 - 4.09        | 47720             | 7360       | 2           | 7351         | 1037         | 30.9        | 0.43          | 0.025           | 0.025           |
| 4.09 - 3.25         | 57550             | 7401       | 4           | 7384         | 947          | 32.2        | 0.47          | 0.029           | 0.027           |
| 3.25 - 2.84         | 31457             | 7430       | 9           | 7409         | 923          | 30.8        | 0.58          | 0.033           | 0.029           |
| 2.84 - 2.58         | 19296             | 7338       | 18          | 7306         | 899          | 29.5        | 0.69          | 0.039           | 0.030           |
| 2.58 - 2.39         | 16529             | 7280       | 14          | 7249         | 884          | 29.1        | 0.76          | 0.040           | 0.031           |
| 2.39 - 2.25         | 14689             | 7178       | 25          | 7129         | 863          | 28.4        | 1.18          | 0.058           | 0.033           |
| 2.25 - 2.14         | 12120             | 7213       | 21          | 7177         | 874          | 27.4        | 1.57          | 0.076           | 0.036           |
| 2.14 - 2.05         | 10195             | 6939       | 37          | 6881         | 844          | 25.7        | 2.05          | 0.092           | 0.038           |
| 2.05 - 1.97         | 7784              | 4574       | 60          | 4453         | 742          | 19.8        | 2.19          | 0.097           | 0.039           |
| 1.97 - 1.90         | 6051              | 1575       | 10          | 1476         | 387          | 14.5        | 2.10          | 0.093           | 0.040           |
| 24.91 - 1.90        | 24681             | 64288      | 200         | 63815        | 8400         | 27.6        | 1.06          | 0.040           | 0.040           |



**Figure 9.** A diffraction image showing the resolution of a 421 Å axis using the MicroMax.

**Table 9.** Processing results for a sample with a 421 Å axis collected with the MicroMax in the long axis configuration.

| Resolution range | Average counts | Num obs | Num rejs | Num ovlp | Num mults | <<I>/<sig>> | ChiSq norm | Rmerge shell | Rmerge cumul |
|------------------|----------------|---------|----------|----------|-----------|-------------|------------|--------------|--------------|
| 42.34 - 7.61     | 2970           | 3954    | 40       | 3528     | 1134      | 20.4        | 1.05       | 0.037        | 0.037        |
| 7.61 - 6.05      | 1257           | 3966    | 16       | 3674     | 1189      | 14.7        | 0.82       | 0.051        | 0.041        |
| 6.05 - 5.28      | 1455           | 4000    | 24       | 3749     | 1223      | 14.8        | 1.03       | 0.061        | 0.046        |
| 5.28 - 4.80      | 2023           | 3973    | 50       | 3712     | 1204      | 16.0        | 1.36       | 0.069        | 0.052        |
| 4.80 - 4.46      | 2494           | 2660    | 53       | 2334     | 699       | 17.3        | 1.57       | 0.070        | 0.055        |
| 4.46 - 4.20      | 2486           | 1745    | 23       | 1455     | 415       | 16.1        | 1.29       | 0.066        | 0.056        |
| 4.20 - 3.99      | 2022           | 1108    | 5        | 860      | 259       | 12.8        | 1.40       | 0.084        | 0.058        |
| 3.99 - 3.81      | 1711           | 727     | 4        | 522      | 188       | 10.2        | 1.27       | 0.081        | 0.058        |
| 3.81 - 3.67      | 1540           | 388     | 1        | 248      | 107       | 7.6         | 1.02       | 0.082        | 0.059        |
| 3.67 - 3.54      | 1387           | 114     | 0        | 20       | 10        | 5.4         | 0.30       | 0.058        | 0.059        |
| 42.34 - 3.54     | 2021           | 22635   | 216      | 20102    | 6428      | 15.9        | 1.16       | 0.059        | 0.059        |